DeskShaper

Deskshaper Basic
BIS515

EN Instructions Manual

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DE Handbuch

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Bisshops BV
Herengracht 338
1016 CG Amsterdam
The Netherlands
www.deskshaper.com

INDEX

EN - Manual	2
NL - Handleiding	3
DE - Handbuch	4
FR - Mode d'emploi	5

Vinz Deskshaper Basic Pedal Exerciser

INSTRUCTIONS

Use the Deskshaper Basic Pedal Exerciser from a desired seating position. Place the DeskShaper Pedal Exerciser on the floor and start cycling. The DeskShaper Basic Pedal Exerciser can be used with both arms and legs. With the black button the resistance of the DeskShaper Basic Pedal Exerciser can be adjusted.

MAINTENANCE

1. keep the Vinz Deskshaper Basic Pedal Exerciser smooth by spraying a little oil time to time.
2. Store the Vinz Deskshaper Basic Pedal Exerciser in a dry place.
3. Clean the Vinz Deskshaper Basic Pedal Exerciser with a mild all-purpose cleaner and dry with a dry cloth.
4. 1 year warranty

WARNING

1. Ensure a firm and level surface, and avoid loose carpets and / or uneven surfaces.
2. For slippery surfaces, use an anti-slip mat so that the Vinz Deskshaper Basic Pedal Exerciser does not move.
3. If applicable, follow the doctor's advice when developing the personal training program.
4. Never use the Vinz Deskshaper Basic Pedal Exerciser with bare feet or socks. Use footwear, such as running shoes or cross training shoes.
5. Sit down at all times before using the Vinz Deskshaper Basic Pedal Exerciser.
6. Pay attention! Allow the Vinz Deskshaper Basic Pedal Exerciser to cool down completely after use before storing or moving.

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